

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 1 of 13

CORRELATIVE SCIENCE PROCEDURE MANUAL

1. Purpose

This document describes the procedures required for the collection, shipping, and processing of biospecimens from all patients enrolled or registered on A022101. This document also describes the procedures that will be followed subsequent to the receipt of biospecimens by the Alliance Biorepository (i.e. Siteman Cancer Center Tissue Procurement Core at Washington University), prior to their use for protocol-specified and future, unspecified correlative science research studies. This document should be used by staff involved with any aspect of the A022101 biospecimen collection, processing, and submission; including staff at satellite institutions.

2. Scope

This document applies to all biospecimens collected specifically for A022101 only. Please refer to the trial protocol-specific language for additional details regarding eligibility, participant enrollment, data submission, and specific procurement procedures. **Please ensure that you are reading the most updated version of this document. This document may experience minor updates, revisions, and clarifications independent of a formal protocol amendment. The most recent version of this document may be found on the Alliance website and CTSU.**

3. Definitions

Term	Definition
ABWUSTL	Alliance Biorepository at Washington University in St. Louis
FFPE	Formalin fixed, paraffin embedded
H&E	Hematoxylin and Eosin

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL	Version No: 1.1	Effective Date: 08/15/2023
	A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Replaces: 1.0	Page 2 of 13

4. Contact Information

Protocol-related questions may be directed as follows:	
Questions	Contact (via email)
Questions regarding patient eligibility, treatment, and dose modification:	Study Chair: Eric D. Miller, MD, PhD eric.miller@osumc.edu Nursing Contact: Barbara Kleiber, RN barbara.kleiber@osumc.edu Protocol Coordinator: Jamie Crawley jcrawley@bsd.uchicago.edu (where applicable) Data Manager: Brandon Bright bright.brandon@mayo.edu
Questions related to data submission, RAVE or patient follow-up:	Data Manager: Brandon Bright bright.brandon@mayo.edu
Questions regarding the protocol document and model informed consent:	Protocol Coordinator: Jamie Crawley jcrawley@bsd.uchicago.edu
Questions related to IRB review:	Alliance Regulatory Inbox regulatory@allianceNCTN.org
Questions regarding CTEP-AERS reporting:	Alliance Pharmacovigilance Inbox pharmacovigilance@alliancencn.org
Questions regarding specimens/specimen submissions:	Alliance Biorepository at Washington University: alliance@email.wustl.edu
Questions regarding drug administration:	Pharmacy Contact: Myounghee Lee, Ph.D. mlee1@umm.edu

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 3 of 13

4.1 For information on using the BiOMS system, please refer to the ‘Help’ links on the BiOMS webpage to access the online user manual, FAQs, and training videos. To report technical problems, such as login issues or application errors, please contact: 1-855-55-BIOMS or bioms@alliancenctn.org. For assistance in using the application or questions or problems related to specific specimen logging, please contact: 1-855-55-BIOMS or bioms@alliancenctn.org.

4.2 For all other questions regarding biospecimen procurement and shipping procedures, please contact the Alliance Biorepository Program Manager: 1-314-747-4402 or alliance@email.wustl.edu.

5. Site Preparation

5.1 Please refer to A022101 protocol document for any specific requirements related to patient enrollment, registration, and regulatory compliance.

5.2 Please ensure that you have appropriate log on credentials and can successfully access the BiOMS application. The BiOMS application is used for logging the collection and shipment of biospecimens to the Alliance Biorepository at Washington University. For training and assistance in using the application or questions or problems related to specific specimen logging, please contact: 1-855-55-BIOMS or bioms@alliancenctn.org.

5.3 Identify a reliable source of dry ice for freezing and shipping biospecimens and a -70 to -90 degree Celsius freezer (“ultralow”) in which frozen biospecimens may be stored prior to shipment.

6. Collection Schema

The following biospecimens are to be collected at each of the time points below. Please refer to individual biospecimen collection and processing methods and specific shipping procedures below.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 4 of 13

Time Point	Biospecimen	Quantity	Collection / Processing Method	Shipping	Notes
For patients who consent to A022101 Biobanking					
Within 14 days after pre-registration	Whole blood for plasma	9 x 1 ml aliquots	Frozen plasma (10.1)	Dry Ice	1, 2, 3
Within 14 days after pre-registration	Whole blood for "buffy coat"	3 aliquots	"Buffy Coat" (10.2)	Dry Ice	1, 2, 3
At Randomization	Fixed tissue block from primary tumor	1 block	Fixed tissue block (9.2)	Ambient	1, 4
At Randomization	H&E stained slide AND Tumor tissue scrolls from primary tumor	1 H&E stained slide AND 10 x 10 micron scrolls	H&E stained slide AND Fixed tissue sections (9.3)	Ambient	1, 4
At Randomization	Fixed tissue block from metastatic site	1 block	Fixed tissue block (9.2)	Ambient	1, 4
At Randomization	H&E stained slide AND Tumor tissue scrolls from metastatic site	1 H&E stained slide AND 10 x 10 micron scrolls	H&E stained slide AND Fixed tissue sections (9.3)	Ambient	1, 4
At Randomization	Whole blood for plasma	9 x 1 ml aliquots	Frozen plasma (10.1)	Dry Ice	1, 3
At Randomization	Whole blood for "buffy coat"	3 aliquots	"Buffy Coat" (10.2)	Dry Ice	1, 3
4 months after randomization (+/- 1 month)	Whole blood for plasma	9 x 1 ml aliquots	Frozen plasma (10.1)	Dry Ice	1, 3
4 months after randomization (+/- 1 month)	Whole blood for "buffy coat"	3 aliquots	"Buffy Coat" (10.2)	Dry Ice	1, 3
8 months after randomization (+/- 1 month)	Whole blood for plasma	9 x 1 ml aliquots	Frozen plasma (10.1)	Dry Ice	1, 3
8 months after randomization (+/- 1 month)	Whole blood for "buffy coat"	3 aliquots	"Buffy Coat" (10.2)	Dry Ice	1, 3

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 5 of 13

1 year after randomization (+/- 1 month)	Whole blood for plasma	9 x 1 ml aliquots	Frozen plasma (10.1)	Dry Ice	1, 3
1 year after randomization (+/- 1 month)	Whole blood for "buffy coat"	3 aliquots	"Buffy Coat" (10.2)	Dry Ice	1, 3
Progression	Whole blood for plasma	9 x 1 ml aliquots	Frozen plasma (10.1)	Dry Ice	1, 3
Progression	Whole blood for "buffy coat"	3 aliquots	"Buffy Coat" (10.2)	Dry Ice	1, 3

Notes:

1. Collection is optional but highly encouraged for all patients. All sites are required to offer biospecimen collection to patients during consent. Please see protocol-specific consent documents.
2. Blood collection at pre-registration is only for patients who enrolled to the study prior to initiating systemic therapy.
3. Peripheral blood (EDTA) 3 x 10 ml to be processed for plasma (9 x 1-1.5 ml aliquots) and "buffy coat," frozen on site and shipped on dry ice.
4. A representative, archived primary tumor tissue block from biopsy or surgery should be submitted, if available. An additional fixed tissue block is requested from site of metastasis, if available. If entire tissue block cannot be submitted, one H&E stained slide **AND** ten (10 um) serial tissue scrolls will be accepted as an alternative. If tissue is limited, please submit H&E and as many tissue scrolls as possible. **BLOCK SUBMISSION IS STRONGLY PREFERRED.**

7. Biospecimen Collection Kits

7.1 Blood Specimens

7.1.1 There are no "kits" provided for submission of blood specimens for this study. Sites are responsible for acquiring materials for collection and shipping of these specimens to the Biorepository.

7.2 Tissue Specimens

7.2.1 There are no "kits" provided for submission of the paraffin block, H&E stained slides, or tissue sections for this study.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 6 of 13

7.2.2 Tissue should be packaged to avoid breakage using a padded envelope or, preferably, a small Styrofoam container.

7.2.3 During warm weather months, fixed tissue should be shipped in an insulated container that contains a refrigerant pack, to avoid heat > 25 degrees C (77 degrees F) that may melt paraffin and damage the tissue specimens.

7.2.4 Please see **Section 11 – Biospecimen Shipping** for specific instructions on shipping to ABWUSTL.

8. Biospecimen Labeling and Tracking

8.1 All research biospecimens (cryovials and tissue bags) MUST be labeled with the Alliance patient ID number, patient initials (Last, First, Middle), the date and time (if applicable) of collection and specimen type (i.e. plasma, “buffy coat”).

8.2 Surgical pathology tissue blocks should not be labeled in any manner. The institutional surgical pathology number (e.g. “S16-1234”) and the individual block identifier (e.g. “A3”) should be readable on the block. Provide a **de-identified copy of the surgical pathology report**, labeled with the Alliance patient ID number, corresponding to the blocks submitted. Please ensure the institutional surgical pathology number and block identifier are maintained on the surgical pathology report. See **section 9** for additional details.

8.3 Label all containers and vials with an indelible, solvent-resistant marker when they are at ambient temperature.

8.4 Do not affix any labels to vials or tubes. Label the collection containers directly with the marking pen.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 7 of 13

8.5 All biospecimens that are collected and sent to the Alliance Biorepository must be **logged and tracked in BioMS**. The BioMS system is a web-based application that tracks the collection and shipping of biospecimens. Once individual biospecimens are logged and ‘shipped’ in the BioMS system, a packing manifest will be created by the system. This manifest must be printed out and must accompany all biospecimen shipments. To become familiar with the BioMS system and for further information about training, access, and use, please contact the BioMS Help desk at: 1-855-55-BIOMS or bioms@alliancencn.org.

8.6 In the event that BioMS cannot be accessed, please complete a BioMS Specimen Log and Shipping Manifest form which can be found here- <http://tinyurl.com/alliance-bioms-contingency>.

9. Tissue Collection

9.1 Overview.

9.1.1 Please refer to protocol-specific instructions for procedures related to actual tissue procurement from individual participants. The method for research tissue procurement (needle core biopsy, sampling of surgically resected tumor) is dependent upon the disease site and the individual patient.

9.1.2 When procuring tissue biospecimens by any method, when possible, avoid tissue that is grossly necrotic, hemorrhagic, fatty, or fibrous. If in doubt, briefly (1 min or less) place the tissue segment in a sterile specimen cup containing physiologic (normal) saline to rinse the tissue. Necrotic, hemorrhagic, and fatty tissue will generally dissolve or float on the surface while tumor and parenchymal tissue will remain intact and sink to the bottom of the cup.

9.2 Diagnostic Pathology Fixed Tissue Blocks

9.2.1 For patients who consent to A022101, one fixed tissue block is requested from the primary tumor. An additional fixed tissue block is requested from site of metastasis, if available. Fixed tissue blocks are requested from biopsy or surgery.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 8 of 13

9.2.2 Any clinical surgical pathology block that is submitted for research studies will not be exhausted or rendered unsuitable for future diagnostic use. Any clinical surgical pathology block that is submitted will be returned within ten working days of written request, when needed for clinical management or clinical trial enrollment for a specific patient. Otherwise, all blocks will be returned to the submitting institution when the trial and correlative science study endpoints have been met.


9.2.3 In the event that an institution will not release fixed tissue blocks, the institution may instead submit 1 H&E stained slide **AND** tissue scrolls from each of the requested blocks as an alternative (see **section 9.3**). **BLOCK SUBMISSION IS STRONGLY PREFERRED.**

9.2.4 During warm weather months, paraffin blocks and slides should be shipped in an insulated container that contains a refrigerant pack, to avoid heat > 25 degrees C (77degrees F) that may melt paraffin and damage the tissue specimens.

9.3 H&E Stained Slide and Fixed Tissue Sections

9.3.1 In cases where an institution is unwilling or unable to submit a tissue block for biobanking, a single H&E stained slide for references and serial tissue sections (scrolls, ribbons, curls) from the same block may be submitted.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 9 of 13



- 9.3.2** Cut and perform routine H&E stain on a single section from the tumor tissue block. See figure above for proper mounting and labeling.
- 9.3.3** Cut a “ribbon” (scroll) of 10 paraffin tissue sections at 10 microns. Place the ribbon of tissue directly into a single microcentrifuge tube or any other suitable container. Do not add any water, PBS, or formalin to the tube. The tube should only contain the cut tissue embedded in paraffin. Do not float the tissue ribbon or sections in a water bath or heater as this may melt the paraffin and ruin the sample. Label the tube of tissue following the guidelines outlined above.
- 9.3.4** During warm weather months, paraffin blocks and slides should be shipped in an insulated container that contains a refrigerant pack, to avoid heat > 25 degrees C (77degrees F) that may melt paraffin and damage the tissue specimens.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 10 of 13

10. Blood Collection Methods

10.1 Plasma Processing

- 10.1.1** Collect 10 ml of whole blood by standard venous phlebotomy technique into each of the purple top (EDTA) tubes. A total of 30 ml of whole blood should be collected into the EDTA tubes (3 x 10 ml). Following collection, invert tubes 10 times.
- 10.1.2** Within 2 hours of collection, spin the vacutainer tubes at room temperature in a clinical centrifuge at 2500 xG for 15 minutes.
- 10.1.3** Carefully remove the plasma layer from each vacutainer tube (~3—5 ml in volume per tube), without touching the white, buffy coat layer, and transfer to new 15 ml conical centrifuge tubes. Keep the vacutainer tubes containing the white, buffy coat layers for white blood cell isolation (**section 10.2**).
- 10.1.4** Spin the centrifuge tubes containing plasma at room temperature in a clinical centrifuge at 2500 xG for 15 minutes.
- 10.1.5** Label 9 cryovials as instructed in **section 8**. Make certain each vial is labeled completely and identically.
- 10.1.6** Carefully remove 9 ml of plasma (without touching the pellet) and divide into nine (9) 2 ml labeled cryovials. Each aliquot should be between 1—1.5 ml in volume.
- 10.1.7** Freeze plasma containing cryovials on dry ice or a -70 to -90 degree Celsius ultralow freezer. Store at -70 to -90 degrees C until ready for shipment on dry ice. Frozen plasma should be shipped to the Biorepository within 30 days of collection. Batch shipment is allowed.

10.2 “Buffy Coat” (White Blood Cell) Processing

- 10.2.1** Follow procedures in **section 10.1** for collecting and processing plasma from EDTA tubes.
- 10.2.2** Label 3 cryovials as instructed in **section 8**.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 11 of 13

10.2.3 After removing the plasma, carefully remove the white, “buffy coat” white blood cell layer, avoiding the red blood cell mass as much as possible.

10.2.4 Transfer the buffy coat layer (approximately 0.2 – 0.5 ml) from EDTA tubes into the labeled cryovials. Immediately freeze the cryovials of buffy coat on dry ice or in liquid nitrogen vapor. Do NOT freeze buffy coat cells by placing a warm tube in a -70 to -90 degree Celsius ultralow freezer. Once completely frozen, the cryovials containing the buffy coat cells may be stored at -70 to -90 degrees C until ready for shipment on dry ice. Frozen buffy coat should be shipped to the Biorepository within 30 days of collection. Batch shipment is allowed.

11. Biospecimen Shipping

11.1 Overview

11.1.1 Frozen plasma and buffy coat aliquots should be placed in a biohazard bag inside of a Styrofoam cooler and covered with 3 to 4 lbs (2 kg) of commercially-prepared dry ice. Pellets or chunks are preferred. Make sure the box is filled with dry ice and the weight of the dry ice is noted on the dry ice label on the outside of the shipping container. It is the local sites’ responsibility to obtain dry ice when shipping frozen specimens. Specimens should be shipped according to IATA guidelines. **Frozen aliquots should be shipped to the Biorepository within 30 days of collection. Batch shipment of frozen aliquots is allowed. If sending specimens from multiple patients within a single shipment, please ensure all specimens are properly labeled and logged in the BioMS system. Specimens from each individual patient should be placed into their own biohazard bag that is clearly labeled with the Alliance patient ID number. The accompanying BioMS manifest should be sealed within each individual bag.**

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 12 of 13

11.1.2 A completed copy of the BioMS packing manifest must accompany all shipments. Do not send specimens without a completed BioMS Packing Manifest or substitute “BioMS Downtime Form.” Biospecimens cannot be accepted without this completed form.

11.1.3 If sending tissue, include a copy of the de-identified surgical pathology report.

11.1.4 **Biospecimens should be shipped Monday—Thursday only. Do not ship on Friday, Saturday, Sunday or the day before a nationally recognized holiday.**

11.2 Shipping to ABWUSTL

11.2.1 Ship container according to IATA guidelines and standard institutional policies via **FedEx priority overnight shipping.**

Ship to:

Alliance Biorepository

c/o Siteman Cancer Center Tissue Procurement Core

Washington Univ. School of Medicine

425 S. Euclid Ave.

Room 5120

St. Louis, MO

63110-1005

Phone: 314-454-7615

12. Biospecimen Receipt and Quality Assurance Measures

12.1 Upon receipt, all biospecimens will be accessioned into the TPC informatics system, OpenSpecimen.

12.2 All biospecimens will be logged, associated, and tracked by the unique patient biopsy control number.

12.3 Each individual biospecimen will receive and be physically labeled with a unique biospecimen identifier, associated with the biopsy control number in the TPC informatics system.

ALLIANCE FOR CLINICAL TRIALS IN ONCOLOGY	CORRELATIVE SCIENCE PROCEDURE MANUAL A Pragmatic Randomized Phase III Trial Evaluating Total Ablative Therapy for Patients with Limited Metastatic Colorectal Cancer: Evaluating Radiation, Ablation, and Surgery [ERASur] Short Title- A022101	Version No: 1.1	Effective Date: 08/15/2023
		Replaces: 1.0	Page 13 of 13

12.4 Upon receipt, all physical biospecimens received will be reconciled with what is recorded on the BioMS packing manifest. Any discrepancies noted will be communicated to the Program Manager who will contact the submitting site for reconciliation.

12.5 Upon receipt, any biospecimen received that is not in appropriate physical condition (broken vials, frozen samples that are thawed, ambient samples that are frozen) will be reported to the Program Manager, who will contact the submitting site for reconciliation.

12.6 Aliquoted biofluids will be stored under liquid nitrogen vapor.

12.7 All biospecimens will remain in storage until additional processing or review is requested in writing by the appropriate protocol PI.

13. Document History

Version	Description and Justification of Change	Author	Effective Date
1.1	Added contact table, changed 'tissue sections' to 'tissue scrolls' for clarity, adde warm weather provisioning for shipping	KL	08/15/2023
1.0	New	AAW	12/12/2022